

*Order*  
**ONE RIPETTE® PRO**  
*And get a*  
**STARTER KIT OF 5**  
**UNITS OF THESE 12**  
**DIFFERENT MULTITIPS**



NO▶	1	2	3	4	5	6	7	8	9	10
0.05 ml	1	1.5	2	2.5	3	3.5	4	4.5	5	5.5
0.10 ml	2	3	4	5	6	7	8	9	10	11
0.20 ml	4	6	8	10	12	14	16	18	20	22
0.50 ml	10	15	20	25	30	35	40	45	50	55
1.00 ml	20	30	40	50	60	70	80	90	100	110
1.25 ml	25	37.5	50	62.5	75	87.5	100	112.5	125	137.5
2.50 ml	50	75	100	125	150	175	200	225	250	275
5.00 ml	100	150	200	250	300	350	400	450	500	550
10.00 ml	200	300	400	500	600	700	800	900	1000	1100
12.50 ml	250	375	500	625	750	875	1000	1125	1250	1375
25.00 ml	500	750	1000	1250	1500	1750	2000	2250	2500	2750
50.00 ml	1000	1500	2000	2500	3000	3500	4000	4500	5000	5500
STEPS▶	48	31	22	18	14	12	10	9	8	7

- Great precision
- High compatibility
- Robust
- Tested quality

40010-0000	Multitips 0,05 ml
40078-0000	Multitips extra 0,1 ml
40078-0001	Multitips extra 0,2 ml
40078-0003	Multitips extra 0,5 ml
40078-0004	Multitips extra 1,0 ml
40010-0002	Multitips 1,25 ml
40078-0005	Multitips extra 2,5 ml
40078-0006	Multitips extra 5,0 ml
40078-0007	Multitips extra 10,0 ml
40078-0008	Multitips extra 12,5 ml
40078-0009	Multitips extra 25,0 ml
40078-0010	Multitips extra 50,0 ml



## Products included:

<b>Cat. No.</b>	<b>Description</b>	<b>Size</b>	<b>Price (EUR)</b>
40078-0000	Multitips extra 0,1 mL	100	68,60
40078-0001	Multitips extra 0,2 mL	100	68,60
40078-0003	Multitips extra 0,5 mL	100	68,60
40078-0004	Multitips extra 1,0 mL	100	68,60
40078-0005	Multitips extra 2,5 mL	100	68,60
40078-0006	Multitips extra 5,0 mL	100	68,60
40078-0007	Multitips extra 10,0 mL	100	68,60
40078-0009	Multitips extra 25,0 mL, w/ adapter	25	19,60
40078-0010	Multitips extra 50,0 mL w/ adapter	25	20,42
40078-0011	Multitips extra 0,1 mL, sterile	100	121,68
40078-0012	Multitips extra 0,2 mL, sterile	100	121,68
40078-0013	Multitips extra 0,5 mL, sterile	100	121,68
40078-0014	Multitips extra 1,0 mL, sterile	100	121,68
40078-0015	Multitips extra 2,5 mL, sterile	100	121,68
40078-0016	Multitips extra 5,0 mL, sterile	100	121,68
40078-0017	Multitips extra 10,0 mL, sterile	100	121,68
40078-0019	Multitips extra 25,0 mL, sterile, w/ adapter	25	32,91
40078-0020	Multitips extra 50,0 mL, sterile, w/ adapter	25	33,73