










labclinics

Pipettor compatibility chart

Pipettor Brand/Model	Biohit™	Biohit Proline™	Biohit Proline Plus™	Brand Transferteppette S™	Brand Transferteppette Electronic™	Capp™	Eppendorf Reference™	Eppendorf Research™	Eppendorf Research Plus™	Finnpipette™	Finnpipette Electronic™	Finnpipette F1™	Gilson Pipetmar™	Gilson Pipetman Ultra™	Hamilton™	Labnet BioPette E™	Nichiryo Nichipet EX™	Nichiryo Oxford Benchmate™	Nichiryo Oxford Multimate™	Socorex Calibra 82.2™	VWR Ergonomic High Performance™	VWR Ultra High Performance™
 LAB10 Series - 10 µl Filter Tip		0.5-10 µl		0.5-10 µl		0.5-10 µl	0.5-10 µl	0.5-10 µl	0.5-10 µl	0.5-10 µl 5-50 µl	1-10 µl	1-10 µl	P2 P10		0.2-2 µl 1-10 µl		0.5-10 µl	0.1-2 µl	0.5-10 µl	1-10 µl		0.1-2 µl 0.5-10 µl
 LAB10XL Series - 10 µl Extended Length Tip		0.5-10 µl		0.5-10 µl		0.5-10 µl	0.1-2.5 µl 0.5-10 µl	0.1-2.5 µl 0.5-10 µl	0.1-2.5 µl 0.5-10 µl	0.5-10 µl 5-50 µl	1-10 µl	1-10 µl	P2 P10		0.2-2 µl 1-10 µl		0.5-10 µl	0.1-2 µl	0.5-10 µl	1-10 µl		0.1-2 µl 0.5-10 µl
 LAB10E Series - 10 µl Eppendorf™ Style Filter Tip		0.5-10 µl		0.5-10 µl		0.5-10 µl	0.1-2.5 µl 0.5-10 µl	0.1-2.5 µl 0.5-10 µl	0.1-2.5 µl 0.5-10 µl	0.5-10 µl 5-50 µl	1-10 µl	1-10 µl	P2 P10		0.2-2 µl 1-10 µl		0.5-10 µl	0.1-2 µl	0.5-10 µl	1-10 µl		0.1-2 µl 0.5-10 µl
 LAB20 Series - 20 µl Filter Tip	M100 M200		10-100 µl 20-200 µl	2-20 µl 10-100 µl 20-200 µl	15-300 µl	5-50 µl 10-100 µl 30-300 µl	2-20, 10-100, 50-200 µl	2-20 µl 10-100 µl 20-200 µl 30-300 µl	2-20 µl 10-100 µl 20-200 µl 30-300 µl	2-20 µl 5-50 µl 20-200 µl 30-300 µl	30-300 µl		P20 P100 P200	U20 U200	2.5-25 µl 10-100 µl 30-300 µl		2-20 µl 10-100 µl 20-200 µl	2-20 µl	5-50 µl 50-300 µl	10-100 µl 20-200 µl	2-20 µl 20-200 µl	2-20 µl 10-100 µl 20-200 µl
 LAB100 Series - 100 µl Filter Tip	M100 M200		20-200 µl	20-200 µl	15-300 µl	5-50 µl 10-100 µl 30-300 µl	50-200 µl	20-200 µl 30-300 µl	50-200 µl	5-50 µl 20-200 µl 30-300 µl	30-300 µl		P200	U200	2.5-25 µl 10-100 µl 30-300 µl		20-200 µl		50-300 µl	20-200 µl	20-200 µl	20-200 µl
 LAB200 Series - 200 µl Filter Tip	M100 M200		20-200 µl	20-200 µl	15-300 µl	5-50 µl 10-100 µl 30-300 µl	50-200 µl	20-200 µl 30-300 µl	50-200 µl	5-50 µl 20-200 µl 30-300 µl	30-300 µl		P200	U200	2.5-25 µl 10-100 µl 30-300 µl		20-200 µl		50-300 µl	20-200 µl	20-200 µl	20-200 µl
 LAB1000 Series - 1000 µl Filter Tip		200-1000 µl		100-1000 µl		100-1000 µl	100-1000 µl	100-1000 µl	100-1000 µl	100-1000 µl 200-1000 µl	100-1200 µl		P1000		100-1000 µl	100-1200 µl	100-1000 µl			100-1000 µl	100-1000 µl	100-1000 µl