


## TECHNICAL DATA SHEET

May, 2014

### PLC Filter Pipette Tips

Labclinics offer also high-quality pipette tips with hydrophobic HDPE filters providing an aerosol barrier and protection against contamination of pipettes and samples in the aspiration process. They are compatible with many common pipette brands. The sterile filter tips cover a range from 0.1 µl to 1 ml and are available in racks.



- > Free of RNase and DNase
- > Minimized fluid retention
- > Compatible with many major pipette brands



Cat. No.	Sample Volume (µl)	Color	Packaging	Sterilized	Pack Size (pcs/pack)	Purchase Unit (pcs)
<b>FTRS10-96 (short)</b>	10	Transparent	Rack	+	96	10 x 96
<b>FTR10-96 (long)</b>	10	Transparent	Rack	+	96	10 x 96
<b>FTR20-96</b>	20	Transparent	Rack	+	96	10 x 96
<b>FTR100-96</b>	100	Transparent	Rack	+	96	10 x 96
<b>FTR200-96</b>	200	Transparent	Rack	+	96	10 x 96
<b>FTR1000-96</b>	1000	Transparent	Rack	+	96	10 x 96