

RAD-CONTROL IRRADIATION TAGS


HOW TO GET STARTED



1 Wipe off any condensation from the surface of the blood bag

2 Apply the tag to a clean, dry area

Tip: Keep part of the tag on the blood bag label for best possible adhesion



3 Press down firmly to effectively adhere the tag

Please note: Once adhered the label takes a few minutes time to fully develop its adhesive strength

HELPFUL TIPS

- Refrain from bending the tag, particularly around the indicator window
- Avoid exposure to direct sunlight and UV sunlight
- Store away from sources of radiation

Rad-Control irradiation tags serve to verify whether irradiation has occurred and should not be considered as irradiation dosimeters.

RAD CONTROL
www.onpointmedical.at